

SERIAL NO.	SHORT FIBER (db)			LONG FIBER (db)			INSTALLED	NOTES
	OUR TEST	THEIR TEST	DIFF	OUR TEST	THEIR TEST	DIFF		
MPS-2300 #001								
MPS-2300 #002	-5.10	-4.50	-0.60	-5.30	-4.50	-0.80	SFY0	
MPS-2300 #003	-4.00	-4.55	0.55	-4.30	-4.45	0.15	SFY1	
MPS-2300 #004	-4.40	-4.50	0.10	-4.90	-4.55	-0.35	SFY2	
MPS-2300 #005	-4.80	-4.55	-0.25	-5.30	-4.65	-0.65	SFY3	
MPS-2300 #006	-4.60	-4.55	-0.05	-4.80	-4.50	-0.30	SFY4	
MPS-2300 #007	-4.80	-4.55	-0.25	-5.00	-4.45	-0.55	SFY5	
MPS-2300 #008	-4.60	-4.45	-0.15	-4.80	-4.45	-0.35	SFY6	
MPS-2300 #009	-5.00	-4.60	-0.40	-5.00	-4.60	-0.40	GLT Lab 6/9/03	
MPS-2300 #010	-4.80	-4.70	-0.10	-4.90	-4.75	-0.15	SFY7	
MPS-2300 #011	-5.00	-4.65	-0.35	-5.20	-4.50	-0.70	SPARE Y	
MPS-2300 #012	-5.60	-4.85	-0.75	-5.70	-4.85	-0.85	SFX00	
MPS-2300 #013	-5.80	-4.60	-1.20	-6.00	-4.65	-1.35	GLT Lab 6/9/03	
MPS-2300 #014	-5.30	-4.80	-0.50	-5.40	-4.65	-0.75	SFX01	
MPS-2300 #015	-5.50	-4.60	-0.90	-5.50	-4.55	-0.95	SFX02	
MPS-2300 #016	-5.40	-4.70	-0.70	-5.40	-4.65	-0.75	SFX03	
MPS-2300 #017	-5.10	-4.80	-0.30	-5.10	-4.85	-0.25	SFX04	
MPS-2300 #018	-6.10	-4.60	-1.50	-6.30	-4.60	-1.70	SFX05	IN-LINE SPLICES, BOTH FIBERS
MPS-2300 #019	-5.80	-4.65	-1.15	-5.90	-4.60	-1.30	SFX06	
MPS-2300 #020	-5.60	-4.45	-1.15	-6.10	-4.75	-1.35	SFX07	
MPS-2300 #021	-6.00	-4.65	-1.35	-5.80	-4.50	-1.30	SFX08	
MPS-2300 #022	-5.10	-4.35	-0.75	-5.90	-4.85	-1.05	SFX09	
MPS-2300 #023	-5.90	-4.65	-1.25	-5.90	-4.75	-1.15	SFX10	
MPS-2300 #024	-5.60	-4.50	-1.10	-6.10	-4.70	-1.40	SFX11	
MPS-2300 #025	-5.40	-4.55	-0.85	-5.60	-4.65	-0.95	SFX12	
MPS-2300 #026	-5.50	-4.70	-0.80	-5.70	-4.75	-0.95	SFX13	
MPS-2300 #027	-5.90	-4.75	-1.15	-6.00	-4.80	-1.20	SFX14	
MPS-2300 #028	-4.80	-4.30	-0.50	-5.80	-4.80	-1.00	SFX15	